


Correction of Dropped Foot using the Odstock Dropped Foot stimulator (ODFS)

Clinical Outline Summary

- Dropped foot is the inability to lift the foot and to advance the leg forward in a stepping motion during walking. Dropped foot is a common impairment that limits walking following damage to the brain or spinal cord (stroke, brain/spinal cord injury & MS)
- Functional Electrical Stimulation (FES) is a technique used to produce contractions in paralysed muscles by the application of small pulses of electrical current to the nerves that supply the muscles. The stimulation is controlled in such a way that the movement produced provides useful function. A Dropped Foot Stimulator is an FES orthosis that produces functional movement in the leg by the application of small pulses of electrical current to the common peroneal nerve (CPN). Stimulation of the CPN improves foot lift and leg excursion during walking.
- Dropped foot stimulators can be used to assist walking in any neurological condition resulting from an upper motor-neuron lesion (damage to the brain or spinal cord).
- Since the introduction of the technology over 20 years ago, patients with stroke, brain injury, multiple sclerosis and other neurological conditions have been shown to benefit significantly from use of the **Odstock Dropped Foot Stimulator (ODFS)** during walking (see attached reference list of published studies) .
- The **Odstock Dropped Foot Stimulator (ODFS)** stimulates the common peroneal nerve using self adhesive skin electrodes positioned on the lower leg. Stimulation is timed to the walking cycle by use of a pressure switch placed in the shoe under the heel and wired to the stimulator. Stimulation begins when the heel is lifted from the ground and ends just after heel strike. Stimulation causes the foot to lift and assists with the excursion of the leg forward in the stepping motion. Stimulation also acts to stabilize the ankle when weight is loaded on the foot on stance phase of gait. See below:



- The ODFS is designed to be a **long term orthosis** and used on a daily basis for walking activities. The device can be donned in the morning and worn all day. Stimulation can be paused at any time when the user is resting.
- Power is provided to the device by a 9 volt battery

Evidence from studies: (please see attached reference list below)

- Increased walking safety due to reduced incidence of tripping
- Increased ankle stability in stance phase of gait
- Increased walking speed
- Reduced effort of walking (less fatigue) & ability to walk longer distances
- Increased confidence and independence when walking due to reduced fear of tripping and falling
- Decreased spasticity in the affected leg & improved range of motion at the ankle

The ODFS Pace can be worn at the waist or the knee



Setting & testing walking parameters in the clinic



Development of the ODFS

- The Odstock Dropped Foot Stimulator (ODFS) was developed at the Department of Medical Physics at Salisbury District Hospital, UK in 1989
- In 2000 ¹⁴Use of the ODFS was recommended as an intervention for neurologically impaired gait in the Royal College of Physicians “Clinical Guidelines on Stroke”.
- To date, in the UK, over 2000 patients have been fitted with an ODFS for dropped foot
- In July 2008, published guidelines by the National Institute of Clinical Excellence (NICE) recommended dropped foot stimulators as an intervention in stroke. {www.nice.org.uk/Guidelines/IPG278}.

- In November 2008, the ODFS received Health Canada approval as a medical device for the treatment of dropped foot and this device is now available to patients in Canada.

References

1. BurrIDGE JH, Taylor PN, Hagan SA, Wood DE, Swain ID. The effect of common peroneal nerve stimulation on quadriceps spasticity in hemiplegia. *Physiotherapy*, 83(2): 82-89, 1997.
2. BurrIDGE J, Taylor P, Hagan S, Wood D, Swain I. (1997) The effects of common peroneal nerve stimulation on the effort and speed of walking: A randomised controlled clinical trial with chronic hemiplegic patients. *Clin Rehabil* 11. 201-210.
3. BurrIDGE JH, Taylor PN, Hagan SA, Swain ID. Experience of the clinical use of the Odstock Drop Foot Stimulator. *Artificial Organs* 21(3): 254-260, 1997
4. Swain ID, Taylor PN, BurrIDGE JH, Hagan SA, Wood DE. Report to the development evaluation committee Common peroneal stimulation for the correction of drop-foot (1996)
5. Taylor P, Mann G, Jolley C, Swain I. Economic Justification for the Odstock Dropped Foot Stimulator (ODFS). ISPO meeting 3rd Nov 2007
6. Taylor PN, BurrIDGE JH, Wood DE, Norton J, Dunkerley A, Singleton C, Swain ID. Clinical use of the Odstock Drop Foot Stimulator - its effect on the speed and effort of walking. *Archives of Physical Medicine and Rehabilitation*, 80: 1577-1583, 1999.
7. Taylor PN, BurrIDGE JH, Wood DE, Norton J, Dunkerley A, Singleton, C, Swain ID. Clinical audit of five years provision of the Odstock Drop Foot Stimulator. *Artificial Organs*, 23(5): 440-442, 1999.
8. Taylor PN. The use of electrical stimulation for correction of dropped foot in subjects with upper motor neuron lesions. *Advances in Clinical Neuroscience and Rehabilitation*, 2(1): 16-18, 2002.
9. Swain ID, Taylor PN. The clinical use of functional electrical stimulation in neurological rehabilitation. *In: Horizons in Medicine 16 – Updates on major clinical advances*. Ed. Franklyn J. Pub. Royal College of Physicians, ISBN 1-86016-233-9, London, pp. 315-322, 2004.
10. Taylor PN, BurrIDGE JH, Wood DE, Norton J, Dunkerley A, Singleton, C, Swain ID. Patient perceptions of the Odstock Drop Foot Stimulator. *Clinical Rehabilitation*, 13: 333-340, 1999.
11. Mann GE, Jolley CJ, Taylor PN. An investigation into the effect of functional electrical stimulation on mobility and quality of life in patients with Multiple Sclerosis. *10th Annual Conference of the International FES Society*, pp. 309-311, Montreal, Canada, July 2005.
12. Esnouf JE, Taylor PN. Does the Canadian Occupational Performance Measure determine if the Odstock Drop Foot Stimulator improves activities of daily living for people with multiple sclerosis? *9th Annual Conference of the International FES Society and 2nd FESnet Conference*, (ISBN 1-85899-191-9), pp. 267-269, Bournemouth, UK, September 2004.
13. Mann GE, Finn SM, Taylor PN. A Pilot study to Investigate the feasibility of electrical stimulation to assist gait in Parkinson's Disease. *Neuromodulation*, 11 (2) 2008 accepted for publication
14. Intercollegiate working party for stroke, (2000) *National clinical guidelines for stroke* London, Royal College of Physicians ISBN 1860 161 200
15. Kottink AI, Hermens HJ, Nene AV, Tenniglo MJ, van der Aa HE, Buschman HP, Ijzerman MJ A randomized controlled trial of an implantable 2-channel peroneal nerve stimulator on walking speed and activity in poststroke hemiplegia. *Arch Phys Med Rehabil*. 2007 Aug;88(8):971-